

Agate Mortar Micro Powder Grinder

❖ Technical Specification Parameters

- A. Diameter of Mortar: 120 mm
- B. Height of Mortar: 34 mm
- C. Speed of Agate Pestle: 20-30 RPM
- D. Speed of Agate Mortar: 10 RPM
- E. Power of Pestle Motor: 25w
- F. Power of Mortar Motor: 10w
- G. Material of Mortar & Pestle: Agate
- H. Grinding Time: 1-30 minutes
- I. L X W X H: 240X240X500 mm
- J. Weight: 120 kg
- K. Voltage: 230 v
- L. Pressure of Pestle: Variable & Manual
- M. Loading of Pestle: Spring Loaded



Ashwini Jain
Ants Ceramics Pvt. Ltd